



Anti-RARA polyclonal antibody (DPABT-H21185)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Rabbit Anti-RARA Polyclonal Antibody
Target	RARA
Immunogen	Peptide
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Antigen Affinity Purified
Conjugate	Unconjugated
Applications	ELISA, WB, IHC
Molecular Weight	462aa,51kd
Cellular Localization	Nucleus
Format	Lyophilized Powder
Size	100 µl, 200 µl
Buffer	PBS with 0.02% Sodium Azide, 0.5%BSA,PH7.2
Preservative	0.02% Sodium Azide
Storage	-20 °C, Avoid freeze / thaw cycles

GENE INFORMATION

Gene Name	RARA retinoic acid receptor, alpha [Homo sapiens]
Official Symbol	RARA
Synonyms	RARA; retinoic acid receptor, alpha; retinoic acid receptor alpha; NR1B1; RAR; RAR-alpha; retinoic acid receptor, alpha polypeptide; nuclear receptor subfamily2group B member 1; retinoic acid nuclear receptor alpha variant 1; retinoic acid nuclear recepto
Entrez Gene ID	5914
Protein Refseq	NP_000955
UniProt ID	P10276
Chromosome Location	17q21.1
Pathway	Acute myeloid leukemia, organism-specific biosystem; Acute myeloid leukemia, conserved biosystem; Adipogenesis, organism-specific biosystem; Gene Expression, organism-specific biosystem; Generic Transcription Pathway, organism-specific biosystem; IL-3 Sig
Function	chromatin DNA binding; enzyme binding; metal ion binding; protein binding; protein domain specific binding; protein heterodimerization activity; protein kinase A binding; protein kinase B binding; receptor activity; receptor binding; retinoic acid binding