



Anti-ZEB1 polyclonal antibody (DPABT-H21190)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Rabbit Anti-ZEB1 Polyclonal Antibody
Target	ZEB1
Immunogen	Peptide
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Purified
Conjugate	Unconjugated
Applications	ELISA, WB
Molecular Weight	1124aa,124kd
Cellular Localization	Nucleus
Format	Lyophilized Powder
Size	200 µl
Buffer	PBS with 0.02% Sodium Azide, PH7.2
Preservative	0.02% Sodium Azide
Storage	-20 °C, Avoid freeze / thaw cycles

GENE INFORMATION

Gene Name	ZEB1 zinc finger E-box binding homeobox2[Homo sapiens]
Official Symbol	ZEB1
Synonyms	ZEB1; zinc finger E-box binding homeobox 1; posterior polymorphous corneal dystrophy4, PPCD3, TCF8, transcription factor 8 (represses interleukin 2 expression); zinc finger E-box-binding homeobox 1; AREB6; BZP; NIL 2 A; ZEB; Zfh1a; negative regul
Entrez Gene ID	6935
Protein Refseq	NP_001121600
UniProt ID	B2RBI8
Chromosome Location	10p11.22
Pathway	Hypertrophy Model, organism-specific biosystem; TGF-beta Receptor Signaling Pathway, organism-specific biosystem; Transcriptional misregulation in cancer, organism-specific biosystem; Transcriptional misregulation in cancer, conserved biosystem;
Function	E-box binding; E-box binding; chromatin binding; double-stranded DNA binding; metal ion binding; protein binding; sequence-specific DNA binding; sequence-specific DNA binding transcription factor activity; transcription coactivator activity; transcription