



Anti-ASH1L polyclonal antibody (DPABT-H21194)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Rabbit Anti-ASH1L Polyclonal Antibody
Target	ASH1L
Immunogen	Peptide
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Antigen Affinity Purified
Conjugate	Unconjugated
Applications	ELISA, WB
Molecular Weight	236aa,25kd
Cellular Localization	Nucleus
Format	Lyophilized Powder
Size	100 µl, 200 µl
Buffer	PBS with 0.02% Sodium Azide, 0.5%BSA,PH7.2
Preservative	0.02% Sodium Azide
Storage	-20 °C, Avoid freeze / thaw cycles

GENE INFORMATION

Gene Name	ASH1L ash1 (absent, small, or homeotic)-like (Drosophila) [Homo sapiens]
Official Symbol	ASH1L
Synonyms	ASH1L; ash1 (absent, small, or homeotic)-like (Drosophila); histone-lysine N-methyltransferase ASH1L; ASH1; ASH1L1; huASH1; KMT2H; ASH1-like protein; lysine N-methyltransferase 2H; absent small and homeotic disks protein2homolog; probable histone-lysine N
Entrez Gene ID	55870
Protein Refseq	NP_060959
UniProt ID	Q9NR48
Chromosome Location	1q22
Pathway	Lysine degradation, organism-specific biosystem; Lysine degradation, conserved biosystem; Tight junction, organism-specific biosystem; Tight junction, conserved biosystem;
Function	DNA binding; histone-lysine N-methyltransferase activity; metal ion binding; methyltransferase activity; transferase activity; zinc ion binding;