



Rabbit anti-HES5 Polyclonal antibody (DPABT-H21195)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Rabbit Anti-HES5 Polyclonal Antibody
Target	HES5
Immunogen	Peptide
Isotype	IgG
Source/Host	Rabbit
Purification	Antigen Affinity Purified
Conjugate	Unconjugated
Applications	ELISA, WB
Molecular Weight	166aa,18kd
Cellular Localization	Nucleus
Format	Lyophilized Powder
Size	100 µl, 200 µl
Buffer	PBS with 0.02% Sodium Azide, 0.5%BSA,PH7.2
Preservative	0.02% Sodium Azide
Storage	-20 °C, Avoid freeze / thaw cycles

GENE INFORMATION

Gene Name	HES5 hairy and enhancer of split 5 (Drosophila) [Homo sapiens]
Official Symbol	HES5
Synonyms	HES5; hairy and enhancer of split 5 (Drosophila); bHLHb38;
Entrez Gene ID	256482
Chromosome Location	1p36