



# Anti-MKI67 polyclonal antibody (DPABT-H21198)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

|                       |                                      |
|-----------------------|--------------------------------------|
| Product Overview      | mouse Anti-MKI67 Polyclonal Antibody |
| Target                | MKI67                                |
| Immunogen             | Peptide                              |
| Isotype               | IgG                                  |
| Source/Host           | Mouse                                |
| Species Reactivity    | Human                                |
| Purification          | Serum                                |
| Conjugate             | Unconjugated                         |
| Applications          | ELISA, WB                            |
| Molecular Weight      | 3256aa,359kd                         |
| Cellular Localization | Nucleus                              |
| Format                | Lyophilized Powder                   |
| Size                  | 200 µl                               |
| Buffer                | Serum with 0.02% Sodium Azide        |
| Preservative          | 0.02% Sodium Azide                   |
| Storage               | -20 °C, Avoid freeze / thaw cycles   |

# GENE INFORMATION

|                     |                                                                                        |
|---------------------|----------------------------------------------------------------------------------------|
| Gene Name           | <a href="#">MKI67 antigen identified by monoclonal antibody Ki-67 [ Homo sapiens ]</a> |
| Official Symbol     | MKI67                                                                                  |
| Synonyms            | MKI67; antigen identified; antigen KI-67; proliferation-related Ki-67 antigen; KIA;    |
| Entrez Gene ID      | <a href="#">4288</a>                                                                   |
| Protein Refseq      | <a href="#">NP_001139438</a>                                                           |
| UniProt ID          | <a href="#">P46013</a>                                                                 |
| Chromosome Location | 10q26.2                                                                                |
| Function            | ATP binding; nucleotide binding; protein C-terminus binding; protein binding;          |