



Anti-NEUROG2 polyclonal antibody (DPABT-H21202)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Rabbit Anti-NEUROG2 Polyclonal Antibody
Target	NEUROG2
Immunogen	Peptide
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Antigen Affinity Purified
Conjugate	Unconjugated
Applications	ELISA, WB, IHC
Molecular Weight	272aa,29kd
Cellular Localization	Nucleus (Potential)
Format	Lyophilized Powder
Size	100 µl, 200 µl
Buffer	PBS with 0.02% Sodium Azide, 0.5%BSA,PH7.2
Preservative	0.02% Sodium Azide
Storage	-20 °C, Avoid freeze / thaw cycles

GENE INFORMATION

Gene Name	NEUROG2 neurogenin 2 [Homo sapiens]
Official Symbol	Neurog2
Synonyms	NEUROG2; neurogenin 2; neurogenin-2; Atoh4; bHLHa8; Math4A; ngn 2; protein atonal homolog 4; class A basic helix-loop-helix protein 8; NGN2; ngn-2; MGC46562;
Entrez Gene ID	63973
Protein Refseq	NP_076924
UniProt ID	Q9H2A3
Chromosome Location	4q25
Function	E-box binding; sequence-specific DNA binding;