



Rabbit anti-DNAJC24 Polyclonal antibody (DPABT-H21211)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Rabbit Anti-DNAJC24 Polyclonal Antibody
Target	DNAJC24
Immunogen	Peptide
Isotype	IgG
Source/Host	Rabbit
Purification	Antigen Affinity Purified
Conjugate	Unconjugated
Applications	ELISA, WB
Molecular Weight	148aa,17kd
Format	Lyophilized Powder
Size	100 µl, 200 µl
Buffer	PBS with 0.02% Sodium Azide, 0.5%BSA,PH7.2
Preservative	0.02% Sodium Azide
Storage	-20 °C, Avoid freeze / thaw cycles

GENE INFORMATION

Gene Name	DNAJC24 DnaJ (Hsp40) homolog, subfamily C, member 24 [Homo sapiens]
Official Symbol	DNAJC24
Synonyms	DNAJC24; DnaJ (Hsp40) homolog, subfamily C, member 24; DPH4, DPH4 homolog (JJJ3, S. cerevisiae) , DPH4, JJJ3 homolog (S. cerevisiae) , ZCSL3, zinc finger, CSL type containing 3; dnaJ homolog subfamily C member 24; JJJ3; 1700030A21Rik; DPH4, JJJ3 homolo
Entrez Gene ID	120526
Protein Refseq	NP_859057
UniProt ID	Q6P3W2
Chromosome Location	11p14.1
Function	heat shock protein binding; metal ion binding; unfolded protein binding;