



Anti-GPR143 polyclonal antibody (DPABT-H21223)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Rabbit Anti-gpr143 Polyclonal Antibody
Target	gpr143
Immunogen	Peptide
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Purified
Conjugate	Unconjugated
Applications	ELISA, WB
Molecular Weight	404aa,44kd
Cellular Localization	Integral membrane protein
Format	Lyophilized Powder
Size	200 µl
Buffer	PBS with 0.02% Sodium Azide, PH7.2
Preservative	0.02% Sodium Azide
Storage	-20 °C, Avoid freeze / thaw cycles

GENE INFORMATION

Gene Name	GPR143 G protein-coupled receptor 143 [Homo sapiens]
Official Symbol	gpr143
Synonyms	GPR143; G protein-coupled receptor 143; OA1, ocular albinism 1 (Nettleship Falls); G-protein coupled receptor 143; ocular albinism 1 (Nettleship Falls); ocular albinism type 1 protein; OA1; NYS6;
Entrez Gene ID	4935
Protein Refseq	NP_000264
UniProt ID	P51810
Chromosome Location	Xp22.3
Pathway	GPCRs, Other, organism-specific biosystem;
Function	G-protein coupled receptor activity; L-DOPA binding; L-DOPA receptor activity; dopamine binding; protein binding; receptor activity; signal transducer activity; tyrosine binding;