



Anti-RPH3AL polyclonal antibody (DPABT-H21236)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Rabbit Anti-RPH3AL Polyclonal Antibody
Target	RPH3AL
Immunogen	Peptide
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Arabidopsis thaliana
Purification	Antigen Affinity Purified
Conjugate	Unconjugated
Applications	ELISA, WB
Molecular Weight	779 aa
Format	Lyophilized Powder
Size	200 µl
Buffer	0.1M Tris-Glycine with 0.02% Sodium Azide, 1%BSA , pH 7.2.
Preservative	0.02% Sodium Azide
Storage	-20 °C, Avoid freeze / thaw cycles

GENE INFORMATION

Gene Name	RPH3AL rabphilin 3A-like (withoutD3domains) [Homo sapiens]
Official Symbol	Rph3al
Synonyms	RPH3AL; rabphilin 3A-like (withoutD3domains); rab effector Noc2; Noc2; noD3domains protein; NOC2;
Entrez Gene ID	9501
Protein Refseq	NP_001177340
UniProt ID	Q9UNE2
Chromosome Location	17p13.3
Function	LIM domain binding; Rab GTPase binding; cytoskeletal protein binding; metal ion binding; zinc ion binding;