



Anti-Dbi polyclonal antibody (DPABT-H21348)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Rabbit Anti-Dbi Polyclonal Antibody
Target	Dbi
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	N/A
Conjugate	Unconjugated
Applications	IB, IHC
Format	Liquid
Buffer	200 µg/ml in PBS with 0.05% NaN ₃ .
Preservative	0.05% Sodium Azide
Storage	Store at 5°C. The antibody can be also aliquotted and stored at -80 °C for long-term storage. Avoid repeated freeze-thawing. Non-hazardous. No MSDS required.

GENE INFORMATION

Gene Name	Dbi diazepam binding inhibitor [Mus musculus]
Official Symbol	Dbi
Synonyms	DBI; diazepam binding inhibitor; acyl-CoA-binding protein; acyl-CoA binding protein; diazepam-binding inhibitor; diazepam binding inhibitor, splice form 1b; EP; Acbp; ACBD1; endozepine;

Entrez Gene ID	13167
Protein Refseq	NP_001033088
UniProt ID	Q4VWZ5
Pathway	PPAR (Peroxisome proliferator-activated receptor) signaling pathway, organism-specific biosystem; PPAR signaling pathway, organism-specific biosystem; PPAR signaling pathway, conserved biosystem;
Function	benzodiazepine receptor binding; binding; fatty-acyl-CoA binding; lipid binding;