



Anti-PROC polyclonal antibody (DPABT-H21456H)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Sheep Anti-Plg PAb
Specificity	Mouse Protein C (blot and ELISA), human Protein C (blot)
Target	Plg
Immunogen	Purified mouse Protein C
Source/Host	Sheep
Species Reactivity	Mouse
Purification	>95% by SDS-PAGE
Conjugate	Unconjugated
Applications	ELISA
Molecular Weight	150000
Format	50% Glycerol/H ₂ O (vol/vol)
Concentration	14.3 mg/mL
Size	1 mg, 5 mg
Buffer	50% Glycerol/H ₂ O (vol/vol)
Preservative	None
Storage	-20 °C

GENE INFORMATION

Gene Name	Plg plasminogen [Mus musculus]
Official Symbol	Plg
Synonyms	PLG; plasminogen; angiostatin; plasmin heavy chain A; plasmin light chain B; Pg; AI649309;
Entrez Gene ID	18815
Protein Refseq	NP_032903
UniProt ID	P20918
Pathway	Blood Clotting Cascade, organism-specific biosystem; Complement and Coagulation Cascades, organism-specific biosystem; Complement and coagulation cascades, organism-specific biosystem; Complement and coagulation cascades, conserved biosystem; Diabetes pathways, organism-specific biosystem; Disease, organism-specific biosystem; Dissolution of Fibrin Clot, organism-specific biosystem;
Function	apolipoprotein binding; catalytic activity; cell surface binding; endopeptidase activity; hydrolase activity; peptidase activity; protein domain specific binding; serine-type endopeptidase activity; serine-type peptidase activity;