



Rabbit anti-Rat Btc Polyclonal antibody (DPABT-H21545H)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Rabbit Anti-Btc PAb
Specificity	As determined by Western Blotting:No cross reactivity with bovine BTC.
Target	Btc
Immunogen	Synthetic epitope of rat beta cellulin (BTC)
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Rat
Conjugate	Unconjugated
Applications	WB
Reconstitution	Dissolve in 500 µl sterile water (a 1/5 dilution). Dilute in a buffer suitable for the desired application
Format	Lyophilized powder.
Size	100 µl
Preservative	None
Storage	At least 2 years at 2-4 °C (lyophilized). After reconstitution store at -20 °C or -70 °C. Avoid freeze-thaw cycles.

GENE INFORMATION

Gene Name	Btc betacellulin [Rattus norvegicus]
Official Symbol	Btc
Synonyms	BTC; betacellulin; probetacellulin;
Entrez Gene ID	64022
Protein Refseq	NP_071592
UniProt ID	Q9JJM4
Pathway	ErbB signaling pathway, organism-specific biosystem; ErbB signaling pathway, organism-specific biosystem; ErbB signaling pathway, conserved biosystem;
Function	epidermal growth factor receptor binding; growth factor activity;