



Anti-CCK polyclonal antibody (DPABT-H21645H)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Rabbit Anti-CCK Polyclonal Antibody
Specificity	CCK-8 100%, CCK-33 100%, CCK-39 84.6%, CCK-8 Nonsulfated
Target	CCK
Immunogen	Synthetic CCK-8 KLH
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	N/A
Conjugate	Unconjugated
Applications	RIA, IHC
Size	50 µl
Preservative	None
Storage	-20 °C, Avoid freeze / thaw cycles

GENE INFORMATION

Gene Name	CCK cholecystokinin [Homo sapiens]
Official Symbol	CCK

Synonyms	CCK; cholecystokinin; MGC117187;
Entrez Gene ID	885
Protein Refseq	NP_000720
UniProt ID	P06307
Chromosome Location	3pter-p21
Pathway	Class A/1 (Rhodopsin-like receptors), organism-specific biosystem; G alpha (q) signalling events, organism-specific biosystem; GPCR downstream signaling, organism-specific biosystem; GPCR ligand binding, organism-specific biosystem; Peptide ligand-binding receptors, organism-specific biosystem; Signal Transduction, organism-specific biosystem; Signaling by GPCR, organism-specific biosystem;
Function	hormone activity; neuropeptide hormone activity;