



Anti-CD79A polyclonal antibody (DPABT-H25779)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Goat Anti-CD79A Polyclonal Antibody
Target	CD79A
Source/Host	Goat
Species Reactivity	N/A
Conjugate	Unconjugated
Size	1 mg
Preservative	None

GENE INFORMATION

Gene Name	CD79A CD79a molecule, immunoglobulin-associated alpha [Canis lupus familiaris]
Official Symbol	CD79A
Synonyms	CD79A; CD79a molecule, immunoglobulin-associated alpha; B-cell antigen receptor complex-associated protein alpha chain; B-cell antigen receptor complex-associated protein alpha-chain; IG-alpha;
Entrez Gene ID	484483
UniProt ID	P0CAN6
Pathway	B Cell Receptor Signaling Pathway, organism-specific biosystem; B cell receptor signaling pathway, organism-specific biosystem; B cell receptor signaling pathway, conserved biosystem;

Primary immunodeficiency, organism-specific biosystem; Primary immunodeficiency, conserved biosystem;
