



Rabbit anti-CAMP Polyclonal antibody (DPABT-H27544)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Rabbit Anti-CAMP Polyclonal Antibody
Target	CAMP
Immunogen	Synthetic peptide corresponding to aa 134-170 of the Human cationic antimicrobial protein 18 (HCAP18).
Isotype	IgG
Source/Host	Rabbit
Purification	Antigen specific affinity Purified
Conjugate	Unconjugated
Size	100 µg
Preservative	None

GENE INFORMATION

Gene Name	CAMP cathelicidin antimicrobial peptide [Homo sapiens]
Official Symbol	CAMP
Synonyms	CAMP; cathelicidin antimicrobial peptide; CAP18; FALL 39; FALL39; LL37; 18 kDa cationic antimicrobial protein; CRAMP; HSD26; CAP-18; FALL-39;
Entrez Gene ID	820

Protein Refseq	NP_004336
UniProt ID	J3KNB4
Chromosome Location	3p21.3
Pathway	Salivary secretion, organism-specific biosystem; Salivary secretion, conserved biosystem; Tuberculosis, organism-specific biosystem; Tuberculosis, conserved biosystem;