



Anti-P2RY12 polyclonal antibody (DPABT-H28327)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Rabbit Anti-P2RY12 Polyclonal Antibody
Target	P2RY12
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	N/A
Conjugate	Unconjugated
Size	50 µg
Preservative	None
Storage	4°C

GENE INFORMATION

Gene Name	P2ry12 purinergic receptor P2Y, G-protein coupled 12 [Mus musculus]
Official Symbol	P2RY12
Synonyms	P2RY12; purinergic receptor P2Y, G-protein coupled 12; P2Y purinoceptor 12; P2Y12; 2900079B22Rik; 4921504D23Rik;
Entrez Gene ID	70839
Protein Refseq	NP_081847

UniProt ID	Q546L4
Pathway	ADP signalling through P2Y purinoceptor 12, organism-specific biosystem; Class A/1 (Rhodopsin-like receptors), organism-specific biosystem; G alpha (i) signalling events, organism-specific biosystem; GPCR downstream signaling, organism-specific biosystem; GPCR ligand binding, organism-specific biosystem; GPCRs, Class A Rhodopsin-like, organism-specific biosystem; Hemostasis, organism-specific biosystem;
Function	ADP receptor activity; ADP receptor activity; G-protein coupled adenosine receptor activity; G-protein coupled purinergic nucleotide receptor activity; G-protein coupled receptor activity; receptor activity; signal transducer activity;