



Anti-GPIHBP1 (C-terminal) polyclonal antibody (DPABT-H30344)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Rabbit Anti-GPIHBP1 Polyclonal Antibody
Specificity	This product is specific for Human HDLBP1/GPIHBP1.
Target	GPIHBP1
Immunogen	This antibody is specific for the C Terminus Region of the target protein.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Affinity purified
Conjugate	Unconjugated
Applications	WB, ELISA, IHC-P
Size	50 µg
Buffer	20 mM Potassium Phosphate, 150 mM Sodium Chloride, pH 7.0
Preservative	None
Storage	Store at -20°C. Avoid freeze-thaw cycles.

GENE INFORMATION

Gene Name	GPIHBP1 glycosylphosphatidylinositol anchored high density lipoprotein binding protein2[Homo sapiens]
Official Symbol	GPIHBP1
Synonyms	GPIHBP1; glycosylphosphatidylinositol anchored high density lipoprotein binding protein 1; GPI anchored high density lipoprotein binding protein 1; glycosylphosphatidylinositol-anchored high density lipoprotein-binding protein 1; GPI HBP1; LOC338328; GPI-
Entrez Gene ID	338328
Protein Refseq	NP_835466
UniProt ID	Q8IV16
Chromosome Location	8q24.3
Function	apolipoprotein binding; chylomicron binding; high-density lipoprotein particle binding; lipase binding; lipid binding; lipoprotein particle binding; protein transmembrane transporter activity;