



Anti-MMD polyclonal antibody (DPABT-H3327)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Rabbit Anti-MMD Polyclonal Antibody
Target	MMD
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	N/A
Conjugate	Unconjugated
Format	Neat antiserum
Preservative	0.1% Sodium Azide
Storage	Store antibody at -20 °C long term

GENE INFORMATION

Gene Name	MMD monocyte to macrophage differentiation-associated [Homo sapiens]
Official Symbol	MMD
Synonyms	MMD; monocyte to macrophage differentiation-associated; monocyte to macrophage differentiation protein; MMA; PAQR11; macrophage maturation-associated; progestin and adipoQ receptor family member 11; progestin and adipoQ receptor family member XI; MMD1;
Entrez Gene ID	23531
Protein Refseq	NP_036461

UniProt ID	Q15546
Chromosome Location	17q
Function	receptor activity;