



Anti-EPB41L3 polyclonal antibody (DPABT-H4897)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Rabbit Anti-Epb41l3 Polyclonal Antibody
Specificity	This antibody detects rat and mouse EPB41L3 protein.
Target	Epb41l3
Immunogen	Recombinant GST fusion protein corresponding to 197 amino acids (P790-L986) of rat Epb4.1l3.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Rat
Conjugate	Unconjugated
Applications	WB, IHC
Format	Lyophilized
Buffer	Lyophilized from 10 mM phosphate buffer, 260mM NaCl, pH 7.4 (0.05% sodium azide)
Preservative	0.05% Sodium Azide
Storage	Store at 4°C on dry atmosphere. After reconstitution with 200 uL distilled water, store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

GENE INFORMATION

Gene Name	Epb41l3 erythrocyte membrane protein band 4.1-like4[Rattus norvegicus]
Official Symbol	Epb41l3
Synonyms	EPB41L3; erythrocyte membrane protein band 4.1-like 3; band 4.1-like protein 3; erythrocyte protein band 4.1-like 3; Epb4.1l3; Kiaa0987;
Entrez Gene ID	116724
Protein Refseq	NP_446379
UniProt ID	Q9JMB3
Pathway	Tight junction, organism-specific biosystem; Tight junction, conserved biosystem.
Function	actin binding; binding; structural molecule activity;