



Anti-CADPS2 (internal region) polyclonal antibody (DPABT-H5434)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Rabbit Anti-CADPS2 Polyclonal Antibody
Target	CADPS2
Immunogen	A synthetic peptide corresponding to internal region 20 amino acids of human CADPS2.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human, Mouse, Rat
Conjugate	Unconjugated
Applications	WB, IHC
Format	Liquid
Size	100 µg
Buffer	In PBS (0.02% sodium azide)
Preservative	0.02% Sodium Azide
Storage	Store at 4°C.

GENE INFORMATION

Gene Name	CADPS2 Ca++-dependent secretion activator 2 [Homo sapiens]
-----------	--

Official Symbol	CADPS2
Synonyms	CADPS2; Ca++-dependent secretion activator 2; Ca++ dependent activator protein for secretion 2; calcium-dependent secretion activator 2; CAPS-2; Ca2+-dependent activator protein for secretion 2; calcium-dependent activator protein for secretion 2; FLJ40851; KIAA1591;
Entrez Gene ID	93664
Protein Refseq	NP_001009571
UniProt ID	Q86UW7
Chromosome Location	7q31.32
Function	lipid binding; metal ion binding;