



Anti-PIWIL4 polyclonal antibody (DPABT-H6215)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Rabbit Anti-PIWIL4 Polyclonal Antibody
Target	PIWIL4
Immunogen	Recombinant protein corresponding to human PIWIL4.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human, Mouse, Rat
Conjugate	Unconjugated
Applications	WB, IHC, ELISA
Format	Liquid
Buffer	In buffer containing 0.02% sodium azide, 2% BSA
Preservative	0.02% Sodium Azide
Storage	Store at -20°C. Aliquot to avoid repeated freezing and thawing.

GENE INFORMATION

Gene Name	PIWIL4 piwi-like 4 (Drosophila) [Homo sapiens]
Official Symbol	PIWIL4

Synonyms	PIWIL4; piwi-like 4 (Drosophila); piwi-like protein 4; FLJ36156; HIWI2; Miwi2; MIWI2; DKFZp686P01248;
Entrez Gene ID	143689
Protein Refseq	NP_689644
UniProt ID	Q7Z3Z4
Chromosome Location	11q12
Pathway	Dorso-ventral axis formation, organism-specific biosystem; Dorso-ventral axis formation, conserved biosystem;
Function	piRNA binding;