



Anti-USP12 (N-terminal) polyclonal antibody (DPABT-H7402)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Rabbit Anti-USP12 Polyclonal Antibody
Target	USP12
Immunogen	A synthetic peptide (conjugated with KLH) corresponding to N-terminus of human USP12.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human, Mouse
Purification	Protein G purification
Conjugate	Unconjugated
Applications	WB, IHC, ELISA
Format	Liquid
Buffer	In PBS (0.09% sodium azide)
Preservative	0.09% Sodium Azide
Storage	Store at 4°C. For long term storage store at -20°C. Aliquot to avoid repeated freezing and thawing.

GENE INFORMATION

Gene Name [USP12 ubiquitin specific peptidase 12 \[Homo sapiens \]](#)

Official Symbol	USP12
Synonyms	USP12; ubiquitin specific peptidase 12; ubiquitin specific protease 12 , ubiquitin specific protease 12 like2, USP12L1; ubiquitin carboxyl-terminal hydrolase 12; ubiquitin thioesterase 12; deubiquitinating enzyme 12; ubiquitin thiolesterase 12; ubiquitin-hydrolyzing enzyme 1; ubiquitin specific protease 12 like 1; ubiquitin-specific-processing protease 12; UBH1; USP12L1;
Entrez Gene ID	219333
Protein Refseq	NP_872294
UniProt ID	O75317
Chromosome Location	13q12.13
Function	cysteine-type endopeptidase activity; peptidase activity; protein binding; ubiquitin thiolesterase activity; ubiquitin thiolesterase activity; ubiquitin-specific protease activity;