



Anti-ASB15 polyclonal antibody (DPABT-H7918)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Rabbit Anti-ASB15 Polyclonal Antibody
Target	ASB15
Immunogen	Recombinant protein corresponding to amino acids of human ASB15.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Antigen affinity purification
Conjugate	Unconjugated
Applications	IHC
Sequence Similarities	YEIVRLLLSHGANVNCYFMHVNDTRFPSVIQYALNDEVMLRLLLNNGYQVEMCFDCMHGD IFGNSFWSEIQEEVLPGWTSCVIKDNPFCEFITVPWMKHLVGRVTRVLIDYMDYVPLCA KLKSALEVQRE
Format	Liquid
Size	100 µl
Buffer	In PBS, pH 7.5 (40% glycerol, 0.02% sodium azide)
Preservative	0.02% Sodium Azide
Storage	Store at 4°C. For long term storage store at -20°C. Aliquot to avoid repeated freezing and

thawing.

GENE INFORMATION

Gene Name	ASB15 ankyrin repeat and SOCS box containing 15 [Homo sapiens]
Official Symbol	ASB15
Synonyms	ASB15; ankyrin repeat and SOCS box containing 15; ankyrin repeat and SOCS box protein 15; FLJ43370; ASB-15; ankyrin repeat and SOCS box-containing 15; DKFZp779M1258;
Entrez Gene ID	142685
Protein Refseq	NP_563616
UniProt ID	Q8WXK1
Chromosome Location	7q31.31
Pathway	Adaptive Immune System, organism-specific biosystem; Antigen processing: Ubiquitination and Proteasome degradation, organism-specific biosystem; Class I MHC mediated antigen processing and presentation, organism-specific biosystem; Immune System, organism-specific biosystem.
