



Anti-C20ORF160 polyclonal antibody (DPABT-H8500)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Rabbit Anti-C20orf160 Polyclonal Antibody
Target	C20orf160
Immunogen	Recombinant protein corresponding to amino acids of human C20orf160.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Antigen affinity purification
Conjugate	Unconjugated
Applications	IHC
Sequence Similarities	IRRLVFPKAGRRAACRSSVSRRLHSMPLYPPDYLDLPQILLCDYLEKEVKFLGHLTWVT SSLNPSSRDELLQLLDTARQLKELPLKTTAEQDSILSLSARCLLLTWDRNEELILRIPTH E
Format	Liquid
Buffer	In PBS, pH 7.5 (40% glycerol, 0.02% sodium azide)
Preservative	0.02% Sodium Azide
Storage	Store at 4°C. For long term storage store at -20°C. Aliquot to avoid repeated freezing and thawing.

GENE INFORMATION

Gene Name	C20orf160 chromosome 20 open reading frame 160 [Homo sapiens]
Official Symbol	C20orf160
Synonyms	C20ORF160; chromosome 20 open reading frame 160; uncharacterized protein C20orf160; dJ310O13.5; FLJ43600;
Entrez Gene ID	140706
Protein Refseq	NP_542192
UniProt ID	Q9NUG4
Chromosome Location	20q11.21