



Anti-DPCR1 polyclonal antibody (DPABT-H9592)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Rabbit Anti-DPCR1 Polyclonal Antibody
Target	DPCR1
Immunogen	Recombinant protein corresponding to amino acids of human DPCR1.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Antigen affinity purification
Conjugate	Unconjugated
Applications	IHC
Sequence Similarities	KTTSTTEKTTRTPEKPTLYSEKTICTKGKNTVPVEKPTENLGNTTLTTETIKAPVKSTEN PEKTA AVTKTIKPSVKVTGDKSLTTTSSHLNKTEVTHQVPTGSFTLITSRTLSSITSEA TGNESHPYLNKDG SQK
Format	Liquid
Buffer	In PBS, pH 7.5 (40% glycerol, 0.02% sodium azide)
Preservative	0.02% Sodium Azide
Storage	Store at 4°C. For long term storage store at -20°C. Aliquot to avoid repeated freezing and thawing.

GENE INFORMATION

Gene Name	DPCR1 diffuse panbronchiolitis critical region2[Homo sapiens]
Official Symbol	DPCR1
Synonyms	DPCR1; diffuse panbronchiolitis critical region 1; diffuse panbronchiolitis critical region protein 1; bCX105N19.6; PBLT; MGC126710; MGC126712;
Entrez Gene ID	135656
Protein Refseq	NP_543146
UniProt ID	E9PEI6
Chromosome Location	6p21.32