



Anti-MMP8 polyclonal antibody [Biotin] (DPABY-525)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	The matrix metalloproteinases (MMPs) consist of 24 known human zinc proteases with essential roles in breaking down components of the extracellular matrix (ECM). Additional MMP substrates include cytokines, chemokines, growth factors and binding proteins, cell/cell adhesion molecules, and other proteinases. With a few exceptions, MMPs share common structural motifs including a pro-peptide domain, a catalytic domain, a hinge region, and a hemopexin-like domain. Synthesized as pro-enzymes, most MMPs are secreted before conversion to their active form. MMP activities are modulated on several levels including transcription, pro-enzyme activation, or by their endogenous inhibitors, tissue inhibitors of metalloproteinases (TIMPs). A subset of MMPs are associated with membranes and designated as membrane-type metalloproteinases (MT-MMP).
Specificity	Detects humanMMP-8 in ELISAs and Western blots. In sandwich immunoassays, less than 0.1% cross-reactivity with recombinant human (rh)rhMMP-1, rhMMP-2, rhMMP-9, rhMMP-13, rhMMP-7, rhMMP-10 and rhMMP-12 is observed.
Immunogen	Mouse myeloma cell line NS0-derived recombinant human MMP-8. Phe21-Gly467 Accession Number AAZ38714
Isotype	IgG
Source/Host	Goat
Species Reactivity	Human
Purification	Antigen Affinity-purified
Conjugate	Biotin
Applications	Western Blot, Immunohistochemistry, ELISA Detection (Matched Pair)

Format	Liquid
Size	50 µg
Buffer	Lyophilized from a 0.2 µm filtered solution in PBS with BSA as a carrier protein.
Preservative	None
Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. 12 months from date of receipt, -20 to -70 °C as supplied. 1 month, 2 to 8 °C under sterile conditions after reconstitution. 6 months, -20 to -70 °C under sterile conditions after reconstitution.

GENE INFORMATION

Gene Name	MMP8 matrix metalloproteinase 8 (neutrophil collagenase) [Homo sapiens (human)]
Official Symbol	MMP8
Synonyms	MMP8; matrix metalloproteinase 8 (neutrophil collagenase); HNC; CLG1; MMP-8; PMNL-CL; neutrophil collagenase; PMNL collagenase; matrix metalloproteinase-8; matrix metalloproteinase 8 (neutrophil collagenase);
Entrez Gene ID	4317
Protein Refseq	NP_002415
UniProt ID	P22894
Chromosome Location	11q22.3
Pathway	Activation of Matrix Metalloproteinases; Collagen degradation; Degradation of the extracellular matrix; Extracellular matrix organization; Matrix Metalloproteinases;
Function	calcium ion binding; metalloendopeptidase activity; serine-type endopeptidase activity; zinc ion binding;