



Goat Anti-Rat ROBO1 polyclonal antibody (DPABY-670-U)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	Detects rat ROBO1 in direct ELISAs and Western blots.
Immunogen	Mouse myeloma cell line NS0-derived recombinant rat ROBO1(aa 19-560)
Isotype	IgG
Source/Host	Goat
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB, IHC Recommended dilution: IHC 5-15 ug/mL WB 0.1 ug/mL
Reconstitution	Reconstitute at 0.2 mg/mL in sterile PBS
Format	Purified, Lyophilized. Reconstitute with 10ul double distilled water.
Buffer	Lyophilized from a 0.2 µm filtered solution in PBS with Trehalose.
Preservative	None
Storage	Store at -20°C.
Ship	Wet ice

GENE INFORMATION

Synonyms	ROBO1; roundabout, axon guidance receptor, homolog 1 (Drosophila); roundabout homolog 1
Entrez Gene ID	58946
UniProt ID	Q55005