



Anti-SOD1 polyclonal antibody (DPABY-730)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	The protein encoded by this gene binds copper and zinc ions and is one of two isozymes responsible for destroying free superoxide radicals in the body. The encoded isozyme is a soluble cytoplasmic protein, acting as a homodimer to convert naturally-occurring but harmful superoxide radicals to molecular oxygen and hydrogen peroxide. The other isozyme is a mitochondrial protein. Mutations in this gene have been implicated as causes of familial amyotrophic lateral sclerosis. Rare transcript variants have been reported for this gene. [provided by RefSeq]
Immunogen	Recombinant fragment corresponding to a region within amino acids 1 and 154 of SOD1 (Uniprot ID#P00441)
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Purified by antigen-affinity chromatography.
Conjugate	Unconjugated
Applications	ELISA Pr*, IHC, WB
Molecular Weight	16 kDa
Cellular Localization	Cytoplasm
Positive Control	A549, HeLa, HepG2, HCT116, mouse brain, Rat brain, *F79, *F123
Format	Liquid
Concentration	0.29 mg/ml

Size	25 µl
Buffer	1XPBS, 1%BSA, 20% Glycerol (pH7). 0.01% Thimerosal was added as a preservative.
Preservative	None
Storage	Keep as concentrated solution. Aliquot and store at -38°C or below. Avoid multiple freeze-thaw cycles.

GENE INFORMATION

Gene Name	SOD1 superoxide dismutase 1, soluble [Homo sapiens (human)]
Official Symbol	SOD1
Synonyms	SOD1; superoxide dismutase 1, soluble; ALS; SOD; ALS1; IPOA; hSod1; HEL-S-44; homodimer; superoxide dismutase [Cu-Zn]; SOD, soluble; indophenoloxidase A; Cu/Zn superoxide dismutase; superoxide dismutase, cytosolic; epididymis secretory protein Li 44;
Entrez Gene ID	6647
Protein Refseq	NP_000445
UniProt ID	P00441
Chromosome Location	21q22.11
Pathway	AGE/RAGE pathway; Amyotrophic lateral sclerosis (ALS); Cellular responses to stress; Detoxification of Reactive Oxygen Species; Dopamine metabolism; FOXA1 transcription factor network; Folate Metabolism; Hemostasis;
Function	Rac GTPase binding; chaperone binding; copper ion binding; identical protein binding; protein binding; protein homodimerization activity; protein phosphatase 2B binding; superoxide dismutase activity; zinc ion binding;