



Anti-Streptomycin polyclonal antibody (DPABY-968)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Streptomycin antibody was raised in sheep using Streptomycin-BSA as the immunogen
Source/Host	Sheep
Species Reactivity	N/A
Conjugate	Unconjugated
Applications	ELISA Pr*
Format	Liquid
Size	100 µl
Buffer	Ig fraction supplied in liquid form in 20 mM phosphate buffer, 150 mM sodium chloride and 0.09% sodium azide, pH 7.2
Preservative	0.09% Sodium Azide
Storage	Store at 2-8 °C for short term storage or -20 °C for up to 2 years

BACKGROUND

Introduction	Streptomycin is a protein synthesis inhibitor. It binds to the small 16S rRNA of the 30S subunit of the bacterial ribosome, interfering with the binding of formyl-methionyl-tRNA to the 30S subunit.
Keywords	Streptomycin; Strep;