



Anti-L. biflexa Polyclonal antibody (DPATB-H81761)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Rabbit Anti-Leptospira biflexa Polyclonal Antibody
Target	L. biflexa
Immunogen	Tissue / cell preparation (Whole organism).
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	L. biflexa
Purification	Protein A purified
Conjugate	Unconjugated
Applications	Conjugation, ICC, IF
Format	Liquid
Size	100 µl
Buffer	0.01M PBS, pH 7.2
Preservative	0.1% Sodium Azide
Storage	Store at +4°C short term. Aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

BACKGROUND

Introduction	Leptospirosis, a contagious disease affecting both animals and humans and spread by infection with a bacterial pathogen called <i>Leptospira</i> , may result in chronic liver and kidney disease and fatality in the dog. <i>Leptospira</i> is a flexible, spiral-shaped, Gra
Keywords	<i>Leptospira biflexa</i> ; Bacteria; Spirochaetes; Spirochaetales; Leptospiraceae; <i>Leptospira</i> ; <i>L. biflexa</i> ; <i>L. australis</i> ; <i>L. canicola</i> ; <i>L. grippityphosa</i> ; <i>L. icterohaemorrhagiae</i> ; <i>L. pomona</i> ; <i>L. serjoe</i>