



# Anti-L. bozemanii Polyclonal antibody (DPATB-H81762)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

<b>Product Overview</b>	Rabbit Anti-Legionella bozemanii Polyclonal Antibody
<b>Target</b>	L. bozemanii
<b>Immunogen</b>	Tissue / cell preparation of Legionella bozemanii (heat inactivated)
<b>Isotype</b>	IgG
<b>Source/Host</b>	Rabbit
<b>Species Reactivity</b>	L. bozemanii
<b>Purification</b>	Whole antiserum
<b>Conjugate</b>	Unconjugated
<b>Applications</b>	Agg
<b>Positive Control</b>	Antigenic suspension.
<b>Format</b>	Liquid
<b>Size</b>	2 ml
<b>Buffer</b>	Whole serum
<b>Preservative</b>	None
<b>Storage</b>	Store at +4°C. Do not freeze.

## BACKGROUND

**Introduction**

Legionella is a Gram-negative bacillus having motility through flagella and is the causative agent of legionellosis. Legionella bozemanii is a bacillus of the genus Legionella that causes pneumonia in humans

---

**Keywords**

Legionella; Legionellaceae; Legionella bozemanii; L. bozemanii

---