



Anti-*L. dumoffii* Polyclonal antibody (DPATB-H81763)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Rabbit Anti- <i>Legionella dumoffii</i> Polyclonal Antibody
Target	<i>L. dumoffii</i>
Immunogen	Tissue / cell preparation of <i>Legionella dumoffii</i> (heat inactivated)
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	<i>L. dumoffii</i>
Purification	Whole antiserum
Conjugate	Unconjugated
Applications	Agg
Positive Control	Antigenic suspension.
Format	Liquid
Size	2 ml
Buffer	Whole serum
Preservative	None
Storage	Store at +4°C. Do not freeze.

BACKGROUND

Introduction	Legionella is a Gram-negative bacillus having motility through flagella and is the causative agent of legionellosis. Legionella dumoffii is implicated in pneumonia.
Keywords	Legionella dumoffii; L. dumoffii; Legionella