



Anti-L. pneumophila type 6 Polyclonal antibody (DPATB-H81768)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Rabbit Anti-Legionella pneumophila serogroup 7 Polyclonal Antibody
Target	L. pneumophila type 6
Immunogen	Tissue / cell preparation of Legionella pneumophila serogroup 7 (heat inactivated).
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	L. pneumophila
Purification	Whole antiserum
Conjugate	Unconjugated
Applications	Agg
Positive Control	Antigenic suspension.
Format	Liquid
Size	2 ml
Buffer	Whole serum
Preservative	None
Storage	Store at +4°C. Do not freeze.

BACKGROUND

Introduction

Legionellapneumophila is a flagellated gram negative bacterium found primarily in warmwater environments. Legionella pneumophila causes a type of pneumonia calledLegionnaire disease and a milder condition called Pontiac fever. Infection isacquired through

Keywords

Legionella pneumophila serogroup 6; Proteobacteria; Gamma Proteobacteria; Legionellales; Legionellaceae; Legionella; L. pneumophila; Legionella pneumophila; Legion
