



# Magic™ Anti-M. synoviae Polyclonal antibody (DPATB-H81797)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

<b>Product Overview</b>	Chicken Anti-Mycoplasma synoviae Polyclonal Antibody
<b>Target</b>	M. synoviae
<b>Immunogen</b>	Mycoplasma synoviae
<b>Isotype</b>	IgY
<b>Source/Host</b>	Chicken
<b>Species Reactivity</b>	M. synoviae
<b>Purification</b>	Whole antiserum
<b>Conjugate</b>	Unconjugated
<b>Applications</b>	ELISA, Neut
<b>Format</b>	Liquid
<b>Size</b>	100 µl
<b>Buffer</b>	Whole serum
<b>Preservative</b>	None
<b>Storage</b>	Store at +4°C short term. Aliquot and store at -20°C or -80°C (add glycerol to a final volume of 50% for extra stability). Avoid repeated freeze / thaw cycles.

## BACKGROUND

**Introduction**

Mycoplasma synoviae is an common pathogen of poultry, with worldwide economical importance, that causes respiratory infections and synovitis. Infection occurs mainly in chickens and turkeys but many other avian species can be affected by this pathogen

---

**Keywords**

M. synoviae; Mycoplasma synoviae; Mycoplasma

---