



Anti-MCM (aa 182-291) polyclonal antibody (DPATB-H81801)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Mouse Anti-MCM Polyclonal Antibody
Target	MCM
Immunogen	Fusion protein: GPFLQSQCCGATENRPVSCPECQTGGPFELNSEQTVYKNY QRVSLQESPGSVPPGRLPRAKDVILLEDLVDSCKPGDEVE LTAIYTHSYESSLNSAHGFP, corresponding to amino acids 182/291 of Macrostromum sp. (marine platyhelminth) MCM. GPFLQSQCCGATENRPVSCPECQTGGPFELNSEQTVYKNY QRVSLQESPG
Isotype	IgG
Source/Host	Mouse
Species Reactivity	Macrostromum
Purification	Whole antiserum
Conjugate	Unconjugated
Applications	WB
Cellular Localization	Nuclear
Format	Liquid
Size	100 µl
Buffer	50% Glycerol
Preservative	None

Storage	Store at +4°C short term. Aliquot and store at -20°C long term. Avoid repeated freeze / thaw cycles.
----------------	--

BACKGROUND

Introduction	The minichromosome maintenance (MCM) 2-7 helicase complex functions to initiate and elongate replication forks.
---------------------	---

Keywords	MCM; minichromosome maintenance (MCM)
-----------------	---------------------------------------
