



Magic™ Anti-Paramyxovirus type 3 Polyclonal antibody (DPATB-H81803)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Chicken Anti-Paramyxovirus Type 3 Polyclonal Antibody
Target	Paramyxovirus type 3
Immunogen	Paramyxovirus Type 3.
Isotype	IgY
Source/Host	Chicken
Species Reactivity	Paramyxovirus
Purification	Whole antiserum
Conjugate	Unconjugated
Applications	ELISA, Neut
Format	Liquid
Size	500 µl
Buffer	Whole serum
Preservative	None
Storage	Store at +4°C short term. Aliquot and store at -20°C or -80°C (add glycerol to a final volume of 50% for extra stability). Avoid repeated freeze / thaw cycles.

BACKGROUND

Introduction

The family of Paramyxoviridae contains viruses that induce a wide range of distinct clinical illnesses in humans. These include measles virus, which in rare instances is followed by subacute sclerosing panencephalitis (SSPE); mumps virus, which has sympto

Keywords

Measles Virus; Measles Virus antigen; Morbillivirus; Paramyxovirus; Paramyxovirus Type 3