



Anti-Methotrexate polyclonal antibody (DPATB-H81809)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Goat Anti-Methotrexate Polyclonal Antibody
Target	Methotrexate
Immunogen	Pure Methotrexate conjugated to hemocyanin
Isotype	IgG
Source/Host	Goat
Species Reactivity	N/A
Purification	Whole antiserum
Conjugate	Unconjugated
Applications	ELISA, RIA
Format	Liquid
Size	100 µl
Buffer	unknown
Preservative	See individual product datasheet
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

BACKGROUND

Introduction

Methotrexate (rINN) abbreviated MTX and formerly known as amethopterin, is an antimetabolite and antifolate drug used in treatment of cancer, autoimmune diseases and as an abortifacient in the induction of medical abortions. It acts by inhibiting the meta

Keywords

METHOTREXATE; METHYLAMINOPTERIN; METHYLAMINOPTERINE; ANTIFOLAN; L-(+)-AMETHOPTERIN; 4-AMINO-N 10-METHYLPTEROYL-L-GLUTAMIC ACID; 4-AMINO-10-METHYLPTEROYLGLUTAMIC ACID; (+)-4-AMINO-10-METHYLFOLIC ACID; 4-AMINO-10-METHYLFOLIC ACID; n-(p-(((2,4-diamino-6-pter