



Anti-Mycobacterium polyclonal antibody (DPATB-H81815)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Rabbit Anti-Mycobacterium Polyclonal Antibody
Target	Mycobacterium
Immunogen	Genus specific antigens extracted from Mycobacterium tuberculosis.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Mycobacterium tuberculosis
Purification	Whole antiserum
Conjugate	Unconjugated
Applications	ELISA
Format	Liquid
Size	200 µl
Buffer	Whole serum
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

BACKGROUND

Introduction

Mycobacterium is a genus of actinobacteria, given its own family, the Mycobacteriaceae. It includes many pathogens known to cause serious diseases in mammals, including tuberculosis and leprosy. Mycobacteria can be classified into several major groups for

Keywords

Mycobacterium; Mycobacteriaceae; symplesiomorphies; cladistics; Bacteria; Actinobacteria; Actinomycetales; Corynebacterineae
