



Magic™ Anti-M. gallisepticum MG Polyclonal antibody (DPATB-H81816)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Chicken Anti-Mycoplasma gallisepticum Polyclonal Antibody
Target	M. gallisepticum MG
Immunogen	Mycoplasma gallisepticum.
Isotype	IgY
Source/Host	Chicken
Species Reactivity	M. gallisepticum
Purification	Whole antiserum
Conjugate	Unconjugated
Applications	ELISA, Neut
Format	Liquid
Size	100 µl
Buffer	Whole serum
Preservative	None
Storage	Store at +4°C short term. Aliquot and store at -20°C or -80°C (add glycerol to a final volume of 50% for extra stability). Avoid repeated freeze / thaw cycles.

BACKGROUND

Introduction

Mycoplasmas are procaryotic microorganisms similar to bacteria but lacking a cell wall. This characteristic makes them more susceptible to environmental factors and, thus, easier to eliminate by the use of disinfectants, heat and sunlight (amongst others)

Keywords

Mycoplasma gallisepticum; MG; Mycoplasma
