



Anti-Myoglobin (native) polyclonal antibody (DPATB-H81818)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Goat Anti-Myoglobin Polyclonal Antibody
Target	Myoglobin
Immunogen	Full length native protein from horse (purified).
Isotype	IgG
Source/Host	Goat
Species Reactivity	Horse
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	ELISA, WB, ICC
Format	Liquid
Size	1 mg
Buffer	PBS, pH 7.2
Preservative	0.1% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

BACKGROUND

Introduction Serves as a reserve supply of oxygen and facilitates the movement of oxygen within muscles.

Keywords MB; MGC13548; MYG_HUMAN; Myoglobin; PVALB;

GENE INFORMATION

Entrez Gene ID [4151](#)

UniProt ID [P02144](#)
