



Anti-Naproxen polyclonal antibody (DPATB-H81824)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Rat Conjugated Naproxen (conjugated) Polyclonal Antibody
Target	Naproxen
Immunogen	Synthetic Naproxen conjugated to protein carriers.
Isotype	IgG
Source/Host	Rat
Species Reactivity	N/A
Purification	Ammonium Sulphate Precipitation
Conjugate	Unconjugated
Applications	ELISA, ICC, IHC-Fr
Format	Liquid
Size	50 µl
Buffer	ddH2O, 50% Glycerol
Preservative	None
Storage	Store at +4°C short term. Aliquot and store at -20°C long term. Avoid repeated freeze / thaw cycles.

BACKGROUND

Introduction	Naproxen sodium (INN) is a nonsteroidal anti-inflammatory drug (NSAID). Naproxen and naproxen sodium are marketed under various trade names, including: Aleve, Anaprox, Antalgin, Apranax, Feminax Ultra, Flanax, Inza, Midol Extended Relief, Nalgesin, Naposi
Keywords	Naproxen; Conjugated Naproxen; (+)-(S)-2-(6-methoxynaphthalen-2-yl)propanoic acid; nonsteroidal anti-inflammatory drug; NSAID