



Anti-Nicarbazin polyclonal antibody (DPATB-H81829)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Sheep Anti-Nicarbazin Polyclonal Antibody
Target	Nicarbazin
Immunogen	Nitrobenzene conjugated to Bovine thyroglobulin
Isotype	IgG
Source/Host	Sheep
Species Reactivity	Bovine
Purification	IgG fraction
Conjugate	Unconjugated
Applications	EIA
Format	Liquid
Size	50 µl
Buffer	20mM PBS, 150mM Sodium chloride, pH 7.2
Preservative	0.09% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

BACKGROUND

Introduction

Nicarbazin has been used in starter rations for several decades as an aid in the prevention of faecal and intestinal coccidiosis in broiler chickens. It may be used in combination with ionophore coccidiostatics. Chemically, it is an equimolar complex of 1

Keywords

Nicarbazin; CL22,791 Urea, N,N"-bis(4-nitrophenyl)-compound with 4,6-dimethyl-2(1H)-pyrimidinone, CycarbN-N"-bis(4-nitrophenyl)urea, compound with 4,6-dimethyl-2(1H)-pyrimidinone (1:1); 4,4"-dinitrocarbanilide ; nicarb; nicoxin; nicarzin