



Anti-Nitrazepam polyclonal antibody (DPATB-H81830)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Sheep Anti-Nitrazepam Polyclonal Antibody
Target	Nitrazepam
Immunogen	Nitrazepam conjugated to BSA
Isotype	IgG
Source/Host	Sheep
Species Reactivity	N/A
Purification	IgG fraction
Conjugate	Unconjugated
Applications	cELISA
Format	Liquid
Size	50 µl
Buffer	20mM PBS, 150mM Sodium chloride, pH 7.2
Preservative	0.09% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

BACKGROUND

Introduction Nitrazepam is a hypnotic drug used in the treatment of moderate to severe insomnia which has sedative and motor impairing properties, as well as anxiolytic, amnestic, anticonvulsant, and skeletal muscle relaxant properties. Nitrazepam is available in 6mg a

Keywords Nitrazepam; Alodorm, Arem, Insoma, Mogadon, Nitrados, Nitrazadon, Ormodon, Paxadorm, Remnos, and Somnite; Pacisyn; 9-nitro-6-phenyl-2,5-diazabicyclo[5.4.0]undeca-5,8,10,12-tetraen-3-one; 7-Nitro-5-phenyl-2,3-dihydro-1H-1,4-benzodiazepin-2-one