



Anti-Nitroserotonin polyclonal antibody (DPATB-H81831)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Rabbit Nitroserotonin (conjugated) Polyclonal Antibody
Target	Nitroserotonin
Immunogen	Chemical / Small Molecule-Nitroserotonin conjugated to a protein carrier.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	N/A
Purification	IgG fraction
Conjugate	Unconjugated
Applications	WB, ELISA, ICC
Format	Liquid
Size	50 µl
Buffer	50% Glycerol, ddH2O
Preservative	None
Storage	Store at +4°C short term. Aliquot and store at -20°C long term. Avoid repeated freeze / thaw cycles.

BACKGROUND

Introduction

Nitroserotonin

Keywords

Nitroserotonin