



Anti-Nitroserotonin polyclonal antibody (DPATB-H81832)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Rabbit Anti-Nitroserotonin Polyclonal Antibody
Target	Nitroserotonin
Immunogen	NO ₂ -Serotonin-Glutaraldehyde-Carriers.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	N/A
Purification	Whole antiserum
Conjugate	Unconjugated
Applications	ELISA
Format	Liquid
Size	50 µl
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

BACKGROUND

Introduction	Nitroserotonin
--------------	----------------

KeywordsNitroserotonin
