



Anti-NO-Glutathione polyclonal antibody (DPATB-H81836)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Rat Anti-NO-Glutathione Polyclonal Antibody
Target	NO-Glutathione
Immunogen	NO-Glutathion-Glutaraldehyde-Carriers.
Isotype	IgG
Source/Host	Rat
Species Reactivity	N/A
Purification	Whole antiserum
Conjugate	Unconjugated
Applications	ELISA
Format	Liquid
Size	50 µl
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

BACKGROUND

Introduction Nitroglutathione is a product of the nitric oxide generating reaction between peroxynitrite and

glutathione.

Keywords

NO-Glutathione
