



Anti-Norescitalopram polyclonal antibody (DPATB-H81843)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Sheep Anti-Norescitalopram Polyclonal Antibody
Target	Norescitalopram
Immunogen	Norescitalopram-BTG
Isotype	IgG
Source/Host	Sheep
Species Reactivity	N/A
Purification	IgG fraction
Conjugate	Unconjugated
Applications	ELISA
Format	Liquid
Size	50 µl
Buffer	20mM PBS, 150mM Sodium chloride, pH 7.2
Preservative	0.09% Sodium Azide
Storage	Shipped at 4°C. Store at -20°C.

BACKGROUND

Introduction	Norescitalopram
Keywords	Norescitalopram