



Anti-Nubbin (internal region) polyclonal antibody (DPATB-H81851)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Rabbit Anti-Nubbin Polyclonal Antibody
Target	Nubbin
Immunogen	Synthetic peptide, derived from an internal region of <i>Publilia modesta</i> Nubbin.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	N/A
Purification	Whole antiserum
Conjugate	Unconjugated
Applications	WB, ICC
Cellular Localization	Nucleus. Predominantly nuclear.
Format	Liquid
Size	200 µl
Buffer	Constituent: Whole serum
Preservative	None
Storage	Store at +4°C short term. Aliquot and store at -20°C long term. Avoid repeated freeze / thaw cycles.

BACKGROUND

Introduction	Specific down-regulator of the NEDD8 conjugation system. Recruits NEDD8, UBD, and their conjugates to the proteasome for degradation. Isoform 1 promotes the degradation of NEDD8 more efficiently than isoform 2.
Keywords	C2H2 and C2HC zinc finger protein; JAGGED LIKE; JGL; NEDD8 ultimate buster 1; Negative regulator of ubiquitin-like proteins 1; NUB; NUB1; NUB1_HUMAN; Renal carcinoma antigen NY-REN-18; T6J4 23; T6J4.23;