



Anti-Obestin polyclonal antibody (DPATB-H81852)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Rabbit Anti-Obestin Polyclonal Antibody
Target	Obestin
Immunogen	Obestin peptide FNAPFDVGILSGVQYQQHSQAL-NH2 conjugated to BSA FNAPFDVGILSGVQYQQHSQAL Run BLAST with Run BLAST with
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Protein G purified
Conjugate	Unconjugated
Applications	RIA, EIA
Format	Liquid
Size	150 µg
Buffer	50% Glycerol, PBS, pH 7.5
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

BACKGROUND

Introduction

Obestin is a brand name of Phentermine, an appetite suppressant, belonging to the amphetamine and phenethylamine class. It stimulates neuron bundles to release the neurotransmitters Catecholamines. The neurotransmitters signal a fight or flight response i

Keywords

Obestin
