



Anti-Ornithine polyclonal antibody (DPATB-H81857)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Rabbit Ornithine (conjugated) Polyclonal Antibody
Target	Ornithine
Immunogen	Synthetic Ornithine conjugated to BSA.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	N/A
Purification	Preadsorbed IgG fraction
Conjugate	Unconjugated
Applications	ELISA
Format	Liquid
Size	50 µl
Buffer	50% Glycerol, ddH2O
Preservative	None
Storage	Store at +4°C short term. Aliquot and store at -20°C long term. Avoid repeated freeze / thaw cycles.

BACKGROUND

Introduction

Ornithine is one of the products of the action of the enzyme arginase on L-arginine, creating urea. Therefore, ornithine is a central part of the urea cycle, which allows for the disposal of excess nitrogen. L-Ornithine is one of the products of the actio

Keywords

Ornithine; L-Ornithine; D-ORNITHINE-OH HCL; D-ORN, HCL; D-ORNITHINE HCL; H-D-ORN-OH HCL; D-ORTHININE MONOHYDROCHLORIDE; D(-)-ORNITHINE HYDROCHLORIDE
